

# KC1 Delta & KC4 Delta

Semi-Automated Coagulation Analysis  
for all Haemostasis clotting assays



## THE FEATURES


- Patented mechanical clot detection for reproducible and reliable results
- Routine and STAT testing using half volumes
- Auto start testing using auto pipette
- KC1 Delta, 1 measuring position with incubation area for 4 samples
- KC4 Delta, 4 measuring positions with incubation area for 8 samples
- Individual cuvettes for all tests

## DATA MANAGEMENT

- Enter calibration curves
- Connect to an optional printer or a PC
- Enter reagent specific information for INR calculations
- Programmable test modes
- Precision monitoring

## APPLICATION

- Full range of clotting assays, PT, aPTT, Fibrinogen, Factors
- Suitable for plasma or whole blood
- Replacement for fibrometer or backup to optical systems
- Maintenance free



PART NUMBER	DESCRIPTION	PACKAGING
G05000	KC1 Delta complete with starter kit (H8 x L21 xW14cm) (1.2kg)	1
N04000	KC4 Delta complete with starter kit (H12 x L45 x W35cm) (6.4kg)	1
<b>Consumables for KC1 Delta &amp; KC4 Delta</b>		
Z04101	Macro univettes for KC1 Delta	220
Z04140	Tetravettes Micro for KC4 Delta	150 X 4
<b>Optional consumables</b>		
Z05100	Macro cuvettes with balls for KC1 Delta	1000
Z01000	Ball Dispenser for Z05100	1
Z05111	Micro cuvettes with balls for KC4 Delta	2000
<b>Optional printers</b>		
Z09160	Printer set KC1 Delta 230v	1
Z09165	Printer set KC4 Delta 230v	1
Z09161	Printer set KC1 Delta and KC4 Delta 110v	1
<b>Printer consumables</b>		
852015	KC Delta Thermal Printer Paper	1

For further information, please contact:



Tcoag Ireland Limited  
 IDA Business Park  
 Bray, Co. Wicklow  
 Ireland  
 Ph.: +353 1 274 32 00  
 Fax: +353 1 274 66 78  
 info@tcoag.com